



Figure 1: Schematic representation of the 1200 bp DNA fragment. The fragment is divided into 12 regions, each with a color-coded box and a label. The labels are: 1. 1000 bp (black), 2. 1000 bp (black), 3. 1000 bp (black), 4. 1000 bp (black), 5. 1000 bp (black), 6. 1000 bp (black), 7. 1000 bp (black), 8. 1000 bp (black), 9. 1000 bp (black), 10. 1000 bp (black), 11. 1000 bp (black), 12. 1000 bp (black). The fragment is also divided into 12 regions, each with a color-coded box and a label. The labels are: 1. 1000 bp (black), 2. 1000 bp (black), 3. 1000 bp (black), 4. 1000 bp (black), 5. 1000 bp (black), 6. 1000 bp (black), 7. 1000 bp (black), 8. 1000 bp (black), 9. 1000 bp (black), 10. 1000 bp (black), 11. 1000 bp (black), 12. 1000 bp (black).



č.v. 07
VÝKRES PŘEDPOKLADANÝCH
ZABORŮ ZPF A PUPEL



SP Studio, s.r.o.,

7 **ИП «Специализированная торговля»**, Белгородская обл., г. Белгород
т./ф. 086719315, e-mail: romashov@spets-torg.ru

	Имя	
	Фамилия, Имя, Отчество	
	Подпись	
	Подпись, печать	

[illegible]
